


Meal Bolus Wizard
 6/1/2016 - 6/14/2016

✕

All Meal Boluses 1


Breakfast - All Boluses

Avg. SG (mg/dL)	SG At Bolus: 143	SG At 2hr: 179
-----------------	------------------	----------------




Lunch - All Boluses

Avg. SG (mg/dL)	SG At Bolus: 138	SG At 2hr: 148
-----------------	------------------	----------------




Dinner - All Boluses

Avg. SG (mg/dL)	SG At Bolus: 168	SG At 2hr: 155
-----------------	------------------	----------------



Night - All Boluses


Avg. SG (mg/dL)	SG At Bolus: 121	SG At 2hr: 233
-----------------	------------------	----------------



Use this report to review your glucose management around mealtimes with your healthcare professional (HCP). Each section for mealtimes are from a timerange you would normally eat these meals. For example, from 6AM to 10AM, all your boluses given during this time frame is accounted for under "Breakfast."

The information above each graph contains average sensor glucose (SG) before (pre-meal) and after (post-meal) a bolus. Remember, not all high glucose events are bad especially after eating, however your glucose values should return to pre-meal range within 3-4 hours after a meal and or bolus.

Bolus Count	1		
Carb Ratio (g/U)	22.0, 18.0		
Avg. Carbs (g)	41 ± 0		
Avg. Bolus (U)	1.9 ± 0.0		



Use this Stats column to understand eating habits.

Are you eating more carbs towards the end of the day? Take a look at your Bolus Count, this is a total number of boluses for the entire reporting period, in this example, it is 14 days. The higher the Bolus Count, the likelihood snacking is happening.

Do you have more than one carb ratio for each meal period? Review them with your HCP for appropriateness.

Do you have boluses in the middle of the night? Does this nighttime snacking affect your morning glucose levels? Discuss these questions with your HCP for possible settings adjustments.




Use this side to take notes and write down questions to your HCP for your next visit. This is a good section to write concerns about mealtime glucose events or symptoms you may experience after taking a bolus. This report can help you participate in your diabetes management with your HCP.

Meal Boluses Analyzed: Boluses with noteworthy pattern or outcome are used for analysis

Time	Analyzed Meals	Observations
AM	Breakfast - Analyzed Meals	
PM	Lunch - Analyzed Meals	
	Dinner - Analyzed Meals	
	Night - Analyzed Meals	

3



Medtronic

18000 Devonshire Street
Northridge, CA 91325
USA
800.646.4633

Toll-free: 800.328.2518
(24-hour technical support for
physicians and medical professionals)

Bolus Wizard, CareLink, MiniMed are registered trademarks of Medtronic MiniMed, Inc.
950M16446-011 20170315 © 2017 Medtronic MiniMed, Inc. All Rights Reserved.

www.medtronicdiabetes.com