

Take the right insulin
dose at the right time[†]

InPen™ smart insulin pen system



Medtronic

Take the guesswork out of insulin dosing

The InPen system has significantly more features compared to a traditional insulin pen, which include:

- Tracks active insulin
- Calculates insulin doses
- Missed dose, long-acting insulin, and blood glucose check reminders
- Insulin expiration & temperature alerts
- Shareable reports for healthcare team
- Real-time[‡] CGM readings with Medtronic CGM
- Integrates with Dexcom[®] G6 and G7 CGM with a three-hour delay
- Connect data from InPen to Apple Health app
- Pair multiple pens to the InPen app
- Travel convenience
- Half-unit delivery



1/2-unit dosing positively affects A1C and blood glucose variability.¹



60% of insulin doses are taken with some insulin on board.²



2 out of 3 people need help calculating their insulin doses.²



A bolus calculator is associated with 0.7-1.0% reduction in A1C.³

InPen Insights Reports


Integrated, reliable diabetes data for easy collaboration with your healthcare provider during your visits.



An insulin pen, but smarter



InPen smart insulin pen is a smart insulin delivery system combining a reusable Bluetooth® insulin pen with an intuitive mobile app to help people with diabetes take the right amount of insulin, at the right time.



Compatible with Novolog®, Humalog®, and Fiasp® cartridges (not included)

Connects via Bluetooth®

Monitors insulin temperature

Battery lasts a full year, with no need to recharge

Delivers half-unit doses

With the InPen system, you can have the freedom of an insulin pen without sacrificing data and decision support. InPen changes

how you manage your multiple daily injections (MDI) by tracking your data and providing personalized dosing recommendations to help you make smarter decisions and achieve better results.¹


Meal therapy modes

InPen takes some of the mental math out of your diabetes management by helping you take the right amount of insulin at the right time.


The InPen app offers three standard modes: Carb counting, meal estimation, or fixed dose.



Carb counting




Meal estimation



Fixed dose

Dose calculator with active insulin tracking

Recommends doses based on blood glucose, carb intake, and active insulin.




Reminders across all devices

Notifications include meal and long-acting dose reminders, bedtime and two-hour post-meal blood glucose check reminders, insulin expiration, and temperature alerts.

The InPen app

Real-time glucose

Calculates personalized doses



Tracks active insulin

Glucose history

Syncs with CGMs

Reminds you to dose

Automatically logs doses given by InPen

Creates shareable reports with up to 90 days' data



Get started today! Scan the QR code,
visit medtronic.com/inpen
or call 800-646-4633

† Compatible smart phone, connectivity, proper settings, and accurate and timely logging of all insulin doses from all sources required. See user guide for complete details.

‡ Data may not appear or be delayed in certain instances, including when there is no internet connection. A confirmatory fingerstick is required prior to making adjustments to diabetes therapy.

Smart insulin pens connect to a mobile app to provide dosing calculations, reminders and CGM system integration.

1. Aanstoot H-J, Rodriguez H, Weinzimer S, Vint N, Koeneman L. Precision Dosing of Rapid-Acting Insulin Matters. *Diabetes Technology & Therapeutics*, 2020;22(5).

2. Zaugg, S.D., et al. Diabetes numeracy and blood glucose control: association with type of diabetes and source of care. *Clinical diabetes: a publication of the American Diabetes Association* 32.4 (2014): 152-157.

3. Anderson DG. Multiple daily injections in young patients using the ezy-BICC bolus insulin calculation card, compared to mixed insulin and CSII. *Pediatric Diabetes*. 2009; 10(5):304-309.

4. Klonoff D, Kerr D. Smart Pens Will Improve Insulin Therapy. *J Dia Sci Tech*, 2018;12(3):551-553.

Humalog® is a registered trademark of Eli Lilly and Company. Fiasp® and Novolog® are registered trademarks of Novo Nordisk A/S. The Bluetooth® Wordmark and logos are owned by Bluetooth® SIG, Inc. App Store® and the App Store logo are trademarks of Apple Inc. Google Play™ and the Google Play logo are trademarks of Google LLC.

Terms and Conditions: Eligibility Requirements for InPen™ Patient Access Program: Open to uninsured patients and patients with commercial prescription insurance and who are not enrolled in any government-funded program that pays for prescription devices or treatments. Excluded participants include patients enrolled in any federal, state, or government-funded healthcare program such as Medicare, Medicare Advantage, Medicare Part D, Medicaid, Medigap, Veterans Affairs (VA), the Department of Defense (DoD) or TRICARE® that provides coverage for InPen, or where prohibited by law. Cash Discount Cards and other non-insurance plans are not valid as primary under this offer. If at any time a patient becomes enrolled under any such federal, state, or government-funded healthcare program, he/she will no longer be able to use this program and must call 844-843-7903 to stop participation. Restrictions may apply. Offer subject to change or discontinuation without notice. This is not health insurance. Offer applies to a maximum of 2 prescribed InPens per Patient per benefit year. Offer good only in the USA at participating pharmacies and cannot be redeemed at government-subsidized clinics.

Important Safety Information: InPen™ system

The InPen™ system consists of an app and home-use reusable pen injector. The system is for single-patient use by people with diabetes under the supervision of an adult caregiver, or by a patient age 7 and older for the self-injection of a desired dose of insulin and for calculating an insulin dose or carbohydrate intake based on user entered data. A healthcare professional must assist in dosage programming of the device prior to use, based on various patient-specific criteria and targets. The InPen™ requires a prescription. For additional product and safety information, please consult the Instructions for Use and bit.ly/InPenSafety.

1 year warranty.

US-DBA-2100057 v4.0 ©2025 Medtronic. Medtronic, MiniMed, and logos are trademarks of Medtronic MiniMed, Inc. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.