

# Medtronic

## MiniMed™ 780G system

### Simplified menu map

From the main menu, any item can be easily selected with as few as **three clicks**.



#### HISTORY & GRAPH

- History
- Sensor Glucose Review
- Graph
- Time in Range



#### SMARTGUARD

- SmartGuard Checklist
- Temp Target
- SmartGuard Settings
- SmartGuard On and Off



#### SOUND & VIBRATION

- Silence Sensor Alerts
- Volume
- Sound
- Vibration
- Alert Settings Shortcut



#### RESERVOIR & SET

- New Reservoir & Set
- New Reservoir Only
- New Set Only
- Fill Cannula



#### INSULIN

- Bolus
- Basal
- Suspend/Resume Basal Delivery
- Delivery Settings Shortcut



#### BLOOD GLUCOSE

- BG



#### STATUS

- Suspend All Delivery
- SmartGuard Checklist
- Pump
- Sensor



#### PAIRED DEVICES

- Pair New Device
- Pair CareLink
- Mobile
- Meter
- Sensor



#### SETTINGS

- Alert Settings
- Delivery Settings
- Device Settings